



**veneer**<sup>®</sup>  
E X T R A

ANTERIORS DESIGNED BY NATURE



*Invented by Dr. Sigal Jacobsen*



*Anteriors by Dr. Jan Hajtó*



## **WELCOME TO UVENEER EXTRA** ANTERIOS DESIGNED BY NATURE

A unique and efficient template kit to create beautiful anterior composite restorations. This kit includes a new variety of sizes and shapes to compliment the original Uvener kit. Each template has a natural surface texture based on the line of natural teeth art models presented by Dr. Jan Hajtó, Germany.

Just press on the composite resin, cure and remove the template. You will be amazed by the results.



**XL**  
(M9) Extra-Large






L  
(F1) Large





**M**  
(F12) Medium





**SQ**  
(F3) Square

Well-suited for teeth with short papillae and pediatric cases when the teeth are not fully erupted.





**UVENEER.COM**  
info@uveneer.com

**ULTRADENT.COM**  
800.552.5512